


Flowtaker

UROFLOWMETRY AT HOME




Innovation for Pelvic Health

INTRODUCING FLOWTAKER®

Flowtaker is a unique, affordable single-patient uroflowmeter for home use. Home uroflowmetry can enhance conventional treatment, giving a better representation of the patient's symptoms in their daily lives¹.

Patients can use Flowtaker in the comfort of their own bathroom, simply by placing the supplied beaker on top of the Flowtaker before each void.

At the end of the recording period, the Flowtaker can be returned by mail or in person to allow the healthcare professional to analyze the results. The Flowtaker software supplied with a box of Flowtakers automatically collects and analyses flow traces creating nomograms, voiding diaries, statistical analysis summary and a customizable report.



IMPROVE ASSESSMENT OF PATIENT'S SYMPTOMS

- Save up to two weeks of uroflows
- Record natural voiding patterns including nocturnal uroflows
- View multiple flow traces with overlay options or nomogram
- Automatically generated voiding diary

INCREASE EFFICIENCY AND REDUCE COSTS

- Reduce missed appointments due to no-shows and failed live uroflows
- Convenient software with patient database, analysis and report options
- Time saving for both your patient and the clinic

PATIENT FRIENDLY

- Voiding in patient's private home environment
- Less traveling and spending time in the waiting room
- Easy to use: no buttons to press, and nothing to write down
- 83% of the patients prefer home uroflowmetry²

¹Neurourol Urodyn. 2007;26(4):474-80 'Flow Measurements ... Make A Right?'

²Eur Urol Suppl 2010;9(2):174, 500 'The accuracy of ...cross over study'

Flowtaker was developed by the Departments of Medical Physics and Urology, Newcastle upon Tyne Hospitals NHS Foundation Trust, supported by a grant from the Wellcome Trust.

The Newcastle upon Tyne Hospitals **NHS**
NHS Foundation Trust



LABORIE International
Tel.: +1 905 612 1170
Email: intlmarketing@laborie.com

LABORIE Europe
Tel.: +31 53 480 3700
Email: info@laborie.com

LABORIE USA
Tel.: 1 800 522 6743
Email: usmarketing@laborie.com

Visit WWW.FLOWTAKER.COM
for more information

