

Simple. Safe. Effective.
Treatment for Urinary Disorders.

Urgent® PC

Neuromodulation System





Simple. Safe. Effective.

The Urgent[®] PC Neuromodulation System treats:

- Urinary urgency
- Urinary frequency
- Urge incontinence

These symptoms are commonly associated with overactive bladder (OAB), interstitial cystitis (IC) and post-radical prostatectomy. Urgent PC may be effective even when previous treatments have failed.




Simple Procedure

- Provides Percutaneous Tibial Nerve Stimulation (PTNS)
- Performed in an office, outpatient or hospital setting with minimal operator training
- Easy to use
 - 30-minute treatment sessions
 - A slim needle electrode, inserted near the tibial nerve, carries electric impulses from a hand-held stimulator to the sacral plexus
 - After the 12 initial treatments, maintenance therapy is continued, using the maximum interval between sessions that maintains symptom relief

Safe Alternative

- Provides similar outcomes to drugs without the associated side-effects¹
 - Increased quality of life
 - Decreased incontinence episodes
 - Decreased urge episodes
- An ideal alternative for patients who
 - Do not wish to take drugs
 - Do not respond to drug therapy
 - Are non-compliant with drug therapy
- Well-tolerated by patients
- Effective without the comorbidities and high cost of surgery
- Ideal for use in combination with drug therapy or other treatments if desired

“...results achieved to date with PTNS are comparable with those seen with anticholinergic pharmacologic therapy.”⁹




Effective Performance

- Patient response rates of up to 80%⁶
- Clinically significant reductions in daytime voiding frequency, nighttime voiding frequency and leakage episodes²⁻⁸
- Results maintained at three year mean follow-up¹⁰
- Effective even when previous treatments have failed³⁻⁸

In a double blind placebo controlled study:

- Greater than 50% reduction in incontinence episodes for 75% of PTNS patients¹¹
- Patients in the placebo group showed no reduction or minimal reduction in their symptoms¹¹
- Patients initially in the placebo group, later treated with PTNS, showed similar results to PTNS treated patients¹¹



“Effective treatment for OAB and NOUR [non-obstructive urinary retention], with stable results at 3-year mean follow-up.”¹⁰



Ordering information

Catalog Number	Product Description
UPC 200	Urgent® PC Stimulator 1 Hand-held Stimulator 1 9-Volt Battery 1 Screwdriver 1 Instructions for Use
UPC 250-12	Urgent PC Lead Sets 12 Lead Sets/box 1 Instructions for Use 1 9-Volt Battery

References:

1. Preyer, O., Gabriel, B., Mailath-Pokorny, M., Doerfler, D., Laml, T., Umek, W., et al. (2007). Peripheral tibial neurostimulation (PTNS) versus Tolterodine in the treatment of women with urge urinary incontinence and urge symptoms. Abstract presentation, IUGA. Mexico.
2. van Balken, M.R., Vandoninck, V., Gisolf K.W.H., Vergunst, H., Kiemeneij, L.A.L.M., Debruyne, F.M.J., et al. (2001). Posterior tibial nerve stimulation as neuromodulative treatment of lower urinary tract dysfunction. *J Urol*, 166, 914-918.
3. Vandoninck, V., van Balken, M.R., Finazzi-Agrò, E., Petta, F., Micali, F., Heesakkers, J.P.F.A., et al. (2003). Percutaneous tibial nerve stimulation in the treatment of overactive bladder: Urodynamic data. *NeuroUrol Urodyn*, 22, 227-232.
4. Klingler, H.C., Pycha, A., Schmidbauer, J., & Marberger, M. (2000). Use of peripheral neuromodulation of the S3 region for treatment of detrusor overactivity: A urodynamic-based study. *Urol*, 56, 766-771.
5. Govier, F.E., Litwiller, S., Nitti, V., Kreder, K.J. Jr., & Rosenblatt, P. (2001). Percutaneous afferent neuromodulation for the refractory overactive bladder results of a multicenter study. *J Urol*, 165, 1193-1198.
6. Stoller, M. (1999). Clinical trials of the SANST™ Afferent Nerve Stimulator to treat urinary incontinence: Results on clinical treatment of more than 90 patients. *Eur Urol*, 35(52), 16.
7. Ruiz, B.C., Outeiriño, X.M.P., Martínez, P.C., Dueñas, E.L., & López, A.L. (2004). Peripheral afferent nerve stimulation for treatment of lower urinary tract irritative symptoms. *Eur Urol*, 45, 65-69.
8. Vandoninck, V., van Balken, M.R., Finazzi-Agrò, E., Petta, F., Caltagirone, C., Heesakkers, J.P.F.A., et al. (2003). Posterior tibial nerve stimulation in the treatment of urge incontinence. *NeuroUrol Urodyn*, 22, 17-23.
9. Cooperberg, M., & Stoller, M. (2005). Percutaneous neuromodulation. *Urol Clin N Am*, 32, 71-78.
10. Cappellano, F., Finazzi-Agrò, E., Giollo, A., Petta, F., Catanzaro, M., Miano, R., et al. (2006). Percutaneous tibial nerve stimulation (PTNS): Results at long term follow-up. Abstract presentation, SIUD National Congress. Italy.
11. Finazzi-Agrò, E., et al. (2006). Percutaneous tibial nerve stimulation (PTNS) in the treatment of urge incontinence: A double blind placebo controlled study. Poster, SIUD National Congress. Italy.



Uroplasty BV

Hofkamp 2
6161 DC Geleen, The Netherlands
Tel: + 31 (0)46 423 79 20, Fax: +31 (0)46 423 79 22
e-mail: info.holland@uroplasty.com

Uroplasty Ltd.

Manchester Business Park, 3000 Aviator Way,
Manchester M22 5TG, UK
Tel: 0800 032 5964, Fax: 0800 032 5965
e-mail: info.uk@uroplasty.com

www.uroplasty.com

© 2007 Uroplasty. All rights reserved.
0525004A 10/07



For complete instructions for use, storage, warnings, indications, contraindications, precautions, adverse reactions and disclaimer of warranties, please refer to the insert accompanying each Urgent PC product.
Urgent and Uroplasty are registered trademarks.
U.S. Patent Number 6,493,588. Additional patents pending.